

CHROMagar™ Malassezia

Click below:

● CHROMagar™ Malassezia

For detection of *Malassezia* spp.


Plate Reading

- *Malassezia furfur*
→ large, pale pink and wrinkled
- Other *Malassezia spp*
including *M. globosa* & *M. restricta*
→ mostly pink to purple

For detection of *Malassezia spp.*

Background

Malassezia is a fungi naturally found on the skin of humans. On normal healthy skin it does not cause infections, but when the environment of the skin is altered, *Malassezia* species are able to cause several cutaneous diseases as severe dermatitis or otitis (inflammation of the skin or ears respectively)

In veterinary field, *Malassezia* is also the cause of various infections detected on animal skins (dermatitis).

1 Medium Performance

EASY AND COLOURFUL DETECTION

Since members of the genus *Malassezia* share similar morphological and biochemical characteristics, the use of traditional culture media for differentiating them based on phenotypic features is not suitable. CHROMagar™ Malassezia was developed with the goal of facilitating not only their detection, but also to improve an algorithm for the differentiation of the most common species (see: Revised Culture-Based System for Identification of *Malassezia* Species, by Takamasa et al. JCM No-2007)

Typical Appearance of the colonies reported in the study:

M. pachydermatis CBS 1879 → large, pale pink & smooth
M. restricta CBS 7877 → small, pink & smooth
M. dermatis JCM11348 and JCM11470 → small, pale pink & smooth
M. slooffiae CBS 7956 → large, pale pink & smooth
M. obtusa CBS 7876 → medium, pink & rough
M. globosa CBS 7966 → small, purple & smooth
M. sympodialis CBS 7222 → large, pale pink & smooth
M. furfur CBS 1878 → large, pale pink & wrinkled

Extracted from: «Revised Culture-Based System for Identification of *Malassezia* Species», by Takamasa et al. JCM No-2007)

Medium Description

| | | |
|-------------------------------------------------------------------------------|---------------------------------------|-------------|
| Powder Base | Total | 56.3 g/L |
| | Agar | 15.0 |
| | Peptones and extracts | 38.0 |
| | Chromogenic mix | 2.8 |
| | Chloramphenicol | 0.5 |
| | Storage at 15/30 °C - pH: 6.3 +/- 0.3 | |
| | Shelf Life | > 18 months |
| Mandatory Supplements to the Powder Base (Not provided by CHROMagar™): | | |
| Glycerol | 2 g | |
| Tween 40 | 10 g | |

| | |
|---------------|----------------------------------------------------------------------|
| Usual Samples | Clinical and veterinary samples (skin, ear canals etc.) |
| Procedure | Direct Streaking. Incubation 72 h at 30-37 °C Aerobic conditions. |

Scientific Publications on this product: available on www.CHROMagar.com
 Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com


Order References

Please use these product references when contacting your local distributor:

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