

ANSR™

For Listeria



產品資訊

敏感度： 1 cfu/analytical unit 102  
 cfu/mL post enrichment  
 檢測時間： 20 min  
 包裝規格： 96 tests / 組  
 套組包含： 96 reaction tubes and caps  
 60 mL lysis buffer  
 96 lysis tubes  
 3 lysis reagent bottles


ANSR™

For Salmonella




產品資訊


敏感度： 1 cfu/analytical unit 104  
 cfu/mL post enrichment  
 檢測時間： 10 min  
 包裝規格： 96 tests / 組  
 套組包含： 96 reaction tubes and caps  
 50 mL lysis buffer  
 96 lysis tubes




Target pathogen DNA



Molecular beacon



Amplified DNA segments



Molecular beacons fluoresce when bound to the amplified segments



ANSR reader with case, power cord.  
Ethernet cable and USB cable



ANSR's alternate plug adapters.



**Simpler. Better. Faster.**